

Medical Plaster

- ▶ Latex-free
- **▶** Practical Smartbox

Features

- Does not stick to the wound
- ▶ Breathable material
- ► Water- and oil-resistant
- ► Permanently elastic
- Latex-free, therefore antiallergic
- ▶ Removable without leaving residue

Benefits of use

- Joints remain freely movable
- Plaster, joint dressing, fingernail and fingertip dressing all-in-one
- ▶ Product adheres only to itself
- No cutting necessary because manually tear off

Areas of application

- Everyday injuries
- Smaller cuts
- Abrasions
- Protection against external abrasive influences to the skin













Medical Plaster

Range of use

- ► Trade
- Household
- ► Kindergarten
- ► Sport
- Hobbies
- Pet Owner

Directions for use

Clean the wound from possible contamination. Apply the **CareLine quick wound dressing** depending on the injury type (joint or tip injury) in accordance with the enclosed instructions. After wrapping the tape is easy to tear and to press firmly on the bandage.

Caution

- In the case of incorrect or improper application, any product warranty becomes void.
- For professional and industrial use only.
- Please observe package insert!

Accessories

Art.-No.

beko working gloves -cut-

905 020 2xx



| Technical data | | | |
|--------------------|--|------------|-----------------|
| Basis | A composite nonwoven fabric | | |
| Waterproof | Blood does not leak and the wound does not become wet. You can wash hands normally or work with wet hands. | | |
| Stretchable | Stops excessive bleeding as a compress. It holds firmly when the joint is moved, does not interfere with work. | | |
| Adhesive-free | Does not stick to the wound or to hair, but only to itself. No adhesive residue on your fingers. The dressing is very popular with children because it can be removed painlessly by rolling. | | |
| Colour | beige/blue | | |
| Storage | not above +25° C, dry | | |
| Contents/Packaging | 2 rolls 25 mm à 4,5 m (in a Smartbox) | 12 pcs/box | 3150 pcs/palett |

This publication supersedes all previous editions and renders the latter invalid.

Currently valid documents (safety data sheet) can be viewed on www.beko-group.com!

Subject to technical changes! Print date: 05/2022 Art.-No. 290 8 9EN

